



## Natural Restorations

## 5 STEP WORKFLOW

*Follow our straightforward steps  
for a seamless experience  
from case submission to shipping*

## STEP 1



### YOU PREP YOUR CASE


For help with your preparations, contact us

## STEP 2

### YOU SEND US YOUR PATIENT DATA

- **SCANS**

- Prep scan + antagonist
- Provisional scan + antagonist




- **PHOTOS**

- Facial smiling picture with temps
- Prep picture with shade tab




- **INFO**

- Desired final shade
- Design specifications such as shape and length



## STEP 3



### WE WORK ON YOUR CASE


## STEP 4

### YOU VALIDATE THE CASE IN TWO MOMENTS

#### DESIGN VALIDATION BEFORE WE MANUFACTURE

Check the following on the html

- Overall design
- Macro texture
- Contact points or interproximal spaces
- Intercuspatation or bite
- Emergence profile or gingival relationship
- Pontic design
- Accuracy with provisional or approved design (copy-paste)



---

#### MANUFACTURING VALIDATION BEFORE WE SHIP

Check the following on the final pictures

- Overall harmony
- Luster/micro texture
- Stain and glaze
- Line angles
- Embrasures
- Color



## STEP 5

WE SHIP YOUR CASE TO YOU

